
Use Cases

for

Lark Detect

Version 1.0 approved

Prepared by

Ang Zhan Phung | Ang Yong Xin | Brenda Ng Xin En
Emmanuelle Vania | Sam Jian Shen | How Mo Xuan


Team CoronaSG

4 February 2020

Table of Contents

Use Case Diagram	iii
Use Case List	iv
Use Case Description	v - xi

1. Use Case Diagram



2. Use Case List

<i>Use Case ID</i>	<i>Primary Actor</i>	<i>Use Cases</i>
1	User	Login
2	User	Register
3	User	Reset Password
4	User	Play
5	Database	Generate Report
6	User, Doctor	View Report

3. Use Case Description

Use Case ID:	1		
Use Case Name:	Login		
Created By:	Brenda Ng	Last Updated By:	Brenda Ng
Date Created:	4/2/20	Date Last Updated:	4/2/20

Actors:	User
Description:	The user logs into the application.
Trigger:	NIL
Preconditions:	<ol style="list-style-type: none"> 1. The user must have a pre-existing account in the database. 2. The user must enter a valid NRIC number and the correct password.
Postconditions:	<ol style="list-style-type: none"> 1. The user must be able to reset their password. 2. The application shows the error message upon unsuccessful login. 3. The application prompts user to register new account if NRIC number does not exist in the database.
Normal Flow:	<ol style="list-style-type: none"> 1. The user opens the application. 2. The application requests user to enter his NRIC number and password. 3. The user enters his NRIC number and password. 4. The database will verify if the user's account information are correct. 5. If the account is verified, the application will display the home page.
Alternative Flows:	<ol style="list-style-type: none"> 1.4: If the account password is invalid <ul style="list-style-type: none"> • Application will display error message "Invalid Password!" • Application returns to the login page. 1.4: If the account NRIC does not exist <ul style="list-style-type: none"> • Application will display error message "Account does not exist. Please enter a correct NRIC number or register an account." • Application returns to the login page.
Exceptions:	NIL
Includes:	NIL
Priority:	7
Frequency of Use:	Very High
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	NIL
Notes and Issues:	NIL

Use Case ID:	2		
Use Case Name:	Register		
Created By:	Emmanuelle Vania	Last Updated By:	Emmanuelle Vania
Date Created:	4/2/20	Date Last Updated:	4/2/20

Actors:	User
Description:	Registration page for first time users
Trigger:	Choosing 'Register' from the Login page
Preconditions:	User must not have a pre-existing database
Postconditions:	A new user record is added to the database
Normal Flow:	<ol style="list-style-type: none"> 1. User enters the page 2. User types in their information such as NRIC, name, date of birth, password. 3. User choose to submit the information 4. System verifies that all required information has been provided by the user and saves them to the database. 5. User is automatically logged in.
Alternative Flows:	AF-S4: If system detects invalid or missing information <ul style="list-style-type: none"> • System prompts user to enter the correct information. • System repeats Step 4.
Exceptions:	NIL
Includes:	NIL
Priority:	7
Frequency of Use:	Low
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	NIL
Notes and Issues:	NIL

Use Case ID:	3		
Use Case Name:	Reset Password		
Created By:	How Mo Xuan	Last Updated By:	How Mo Xuan
Date Created:	4 February 2020	Date Last Updated:	4 February 2020

Actors:	User
Description:	User can reset the password if they forget the Lark Detect login's password
Trigger:	
Preconditions:	
Postconditions:	1. User able to login into the game with his/her new password
Normal Flow:	<ol style="list-style-type: none"> 1. User clicks the 'Forget password' button. 2. User received a change password by verifying their personnel information 3. User input the new password in 'New password' and 'Re-enter new password' columns. 4. Database update the password.
Alternative Flows:	NIL
Exceptions:	NIL
Includes:	NIL
Priority:	3
Frequency of Use:	Low
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	NIL
Notes and Issues:	NIL

Use Case ID:	4		
Use Case Name:	Play		
Created By:	Sam Jian Shen	Last Updated By:	Sam Jian Shen
Date Created:	4 February 2020	Date Last Updated:	4 February 2020

Actors:	User
Description:	User must be able to successfully finish the Trail Making Test (TMT) which consist of 2 different tests; test A and test B
Trigger:	After user click “play icon” in main page,
Preconditions:	<ol style="list-style-type: none"> 1. User must successfully login into the system 2. User’s previous attempt 3. User’s profile must specify a valid age. 4. User’s must not failed more than 3 times
Postconditions:	<ol style="list-style-type: none"> 1. User must be able to view the result from the test 2. User result from the test must be saved into the external database system.
Normal Flow:	<ol style="list-style-type: none"> 4.1 The “LarkDetect” system must display the honesty agreement instruction to the user. 4.2 User’s clicked “Agreed” icon. 4.3 The “LarkDetect” system must display the 1st test instruction to the user. 4.4 User’s clicked “Next” icon. 4.5 The “LarkDetect” system must populate with the environment (background/sound effect/assets) and elements (shapes/text/size/font) for the trail of the 1st test. 4.6 The “LarkDetect” system must display a trail of the 1st test to UI. 4.7 The “LarkDetect” system must display the link between the correct sequences of nodes in green colour. 4.8 User’s clicked “Start Game” icon. 4.9 The “LarkDetect” system must populate with the environment (background/sound effect/assets) and elements (shapes/text/size/font) for the 1st test. 4.10 The “LarkDetect” system must display the 1st test to UI 4.11 The “LarkDetect” system must display the link between the correct sequences of nodes in green colour. 4.12 The system must check if 25 sequences of nodes is complete. 4.13 The “LarkDetect” system must display the 2nd test instruction to the user. 4.14 User’s clicked “Next” icon. 4.15 The “LarkDetect” system must populate with the environment (background/sound effect/assets) and elements (shapes/text/size/font) for the trail of the 2nd test. 4.16 The “LarkDetect” system must display a trail of the 2nd test to UI. 4.17 The “LarkDetect” system must display the link between the correct sequences of nodes in green colour. 4.18 User’s clicked “Start Game” icon.

	<p>4.19 The “LarkDetect” system must populate with the environment (background/sound effect/assets) and elements (shapes/text/size/font) for the 2nd test.</p> <p>4.20 The “LarkDetect” system must display the 2nd test to UI</p> <p>4.21 The “LarkDetect” system must display the link between the correct sequences of nodes in green colour.</p> <p>4.22 The system must check if 25 sequences of nodes is complete.</p> <p>4.23 The “LarkDetect” system will store the user’s “performance” to the external system</p> <p>4.24 The “LarkDetect” system will display the result to the user.</p>
Alternative Flows:	<p>4.1 If User’s click the wrong sequence of nodes.</p> <p>4.1.1 The system must display the sequence of link in red color.</p> <p>4.2 If User’s click “disagreed” with honesty agreement.</p> <p>4.2.1 The “LarkDetect” system will pop-up a message-box stating for confirmation “Are you sure?”</p> <p>4.2.2 Move to UC-01</p> <p>4.3 If User’s time taken for the test is more than 3 minutes and failed more than or equal to 3 tries.</p> <p>4.3.1 The system will sign-out the user</p> <p>4.3.2 Move to UC-01</p> <p>4.4 If User’s time taken for the test is more than 3 minutes and failed less than 3 tries</p> <p>4.4.1 The system will move to main page step 4.3 of normal flow</p>
Exceptions:	NIL
Includes:	1
Priority:	10
Frequency of Use:	High
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	NIL
Notes and Issues:	NIL

Use Case ID:	5		
Use Case Name:	Generate Report		
Created By:	Edmund Ang	Last Updated By:	Edmund Ang
Date Created:	4/2/20	Date Last Updated:	4/2/20

Actors:	Database
Description:	To fetch the data from database to generate report
Trigger:	
Preconditions:	1. User will need to play before generating the report
Postconditions:	2. Report and information will be recorded in the database
Normal Flow:	User 5.0 User click on generate the report 5.1 System display the report
Alternative Flows:	NIL
Exceptions:	NIL
Includes:	
Priority:	7
Frequency of Use:	Often
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	Completion of Play case
Notes and Issues:	NIL

Use Case ID:	6		
Use Case Name:	View Report		
Created By:	Ang Yong Xin	Last Updated By:	Ang Yong Xin
Date Created:	4/2/20	Date Last Updated:	4/2/20

Actors:	Doctors Users
Description:	This use case allows both user and doctors to view their history record to make diagnosis on the dementia condition
Trigger:	
Preconditions:	NIL
Postconditions:	Users/Doctors able to view the history report
Normal Flow:	Users: 3. User click on "View Report" button 4. System retrieve user history report from database 5. System displays the history report. Doctors: 6. Doctor click on "View Patient Report" button 7. Doctor enter the patient name into the system 8. System retrieve patient history reports from the database System displays patient history reports.
Alternative Flows:	6.2: No user history 9. System display "No game history" message. 10. Use case return to step 1 5.3: No patient history 11. System display "No game history" message Use case return to step 1
Exceptions:	
Includes:	
Priority:	5
Frequency of Use:	Often
Business Rules:	NIL
Special Requirements:	NIL
Assumptions:	NIL
Notes and Issues:	NIL